

SPECIFICATIONS

ARTICLE

Designation : **Pure Soft disinfectant cleaner**

Code bobet : **24632 : 750ml**

Selling Unit : **Unit**

Non-contractual pictures

**CHARACTERISTICS**

- Action: Fragrance-free disinfectant cleaner
- Appearance : clear, colourless liquid
- pH: 10.75 +/- 0.75
- Density: 1g/cm³ +/- 0.0
- 750 ml can
- Non-corrosive, non-harmful

DESCRIPTION :

- Cleans and disinfects all surfaces
- Formulated without phosphates, allergens, CMR and chlorine

Storage: For optimum shelf life up to the expiry date, keep the product away from heat, light and frost.

Active ingredient: Didecylmethylammonium chloride (CAS N°: 7172-51-5) at 4.5g/kg.

Contains among other ingredients (regulation (EC) N°648/2004): less than 5% non-ionic surfactants; less than 5% EDTA and salts; Disinfectants and didecylmimonium chloride.

INSTRUCTIONS FOR USE :

- Spray evenly over the surface to be cleaned.
- Leave to work for the desired disinfectant effect (see table), then wipe with a clean cloth.

PACKAGING :

Manual sprayers: box of 6 x 750mL

STANDARDS :

Complies with legislation relating to cleaning products for equipment that may come into contact with foodstuffs (Order of 08/09/1999 and its amendments).

TP4: Disinfectants for surfaces in contact with food and animal feed.

TP2: Disinfectants used in the private and public health sectors and other biocides.

NORMES *en condition de saleté générale	SOUCHES	TEMPS DE CONTACT (20°C)
Bactéricide selon les normes EN 1276	<i>Pseudomonas aeruginosa</i> - <i>Escherichia coli</i> - <i>Staphylococcus aureus</i> – <i>Enterococcus hirae</i> – <i>Salmonella enterica</i> - <i>Listeria monocytogenes</i>	5 min
Bactéricide selon les normes EN 13697		15 min
Levuricide- fongicide selon les normes EN1650 et EN13697	<i>Candida Albicans</i> – <i>Aspergillus brasiliensis</i>	15 min
Virucide selon la norme EN14476	<i>Adenovirus type 5</i> - <i>Murine norovirus</i> – <i>Poliovirus</i>	15 min
Virucide selon la norme EN14476	<i>Coronavirus</i>	5 min
Virucide selon la norme EN16777	<i>Adenovirus type 5</i> - <i>Murine norovirus</i> – <i>Poliovirus</i>	15 min