


SPECIFICATIONS



ARTICLE

<u>Designation</u>: Purple scouring pad

Code bobet : 15787 Selling Unit : Unit



Non-contractual pictures

CHARACTERISTICS

"Special Washing Up" pad for intensive scouring, without scratching, replaces stainless steel scouring balls and steel wool

- Composition : synthetic fiber, abrasive grains, synthetic resin

Shape : exclusive butterflyAbrasion : 0g removed

Colour : purpleThickness : 13 mm

- Dimensions: 102 x 133 mm

- Weight: 10.5 g

DESCRIPTION:

This purple butterfly-shaped pad for cleaning and multi-purpose scouring, in high quality synthetic fibers is made of a non woven patch with an open structure impregnated with a resistant resin with abrasive particles homogeneously spread and weakly scratching.

BENEFITS

- Excellent scouring properties, loosen burnt food residues
- Unique butterfly shape, flexible, enables scouring in corners and hard to reach places
- Patch with open structure, easy to clean
- Made in a resin resistant to hot water and detergents and current cleaning products
- Very good durability
- Steady performance throughout the lifetime of the pad
- Perfect substitute to steel wool and stainless steel balls which are identified as "critical points" as they can leave metal particles dangerous if ingested and which are breeding grounds for bacteria

APPLICATION / MAINTENANCE :

- This product can be used wet or dry to scour dishes and pans, clean roasters and most kitchen utensils.
- On steel, stainless steel, kitchen utensils, hotplates, porcelain, glass, ceramic.

IMPORTANT:

- Always carry out tests on an unexposed surface, especially on fragile surfaces.
- Do not use on anti-stick surfaces and treatments, plastics, plexiglass, varnish.
- Do not use on hot surfaces (dry use, maximum temperature 120°C)



MAINTENANCE

- Rinse the patch with clear water after each use.
- If required, remove residues inlaid in the pad with a nylon brush or dip it in a cleaning solution to dissolve burnt food residues.
- Can be more easily washed in a washing machine.
- For perfect disinfection, can be stored in a solution weakly chlorinated (maximum concentration <2.6°ch)
- For better hygiene, renew the pad after 2 or 3 weeks of use.

PACKAGING:

Bag of 24 pads.